## **ONETOUCH**

Select<sup>\*\*</sup>

Blood Glucose Monitoring System Sistema de control de glucosa en sangre

Easy to Use. Results You Can Trust.

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### Owner's Booklet Manual del propietario

#### Welcome to the OneTouch® family!

 $\mathsf{OneTouch}^{\$}$  is committed to creating a world without limits for people with diabetes.

The OneTouch<sup>®</sup> Select<sup>™</sup> Meter makes it simple to get a test result and review past results. This owner's booklet will help you learn how to use your meter properly. Please read it carefully.

If you need assistance, please contact LifeScan Customer Service at 1 800 227-8862 (available 24 hours a day, seven days a week). Please also visit www.LifeScan.com for more information.

**NOTE:** Before testing, read this owner's booklet carefully.

#### **Table of Contents**

Setting up your system	Pg. 1
Starting the test process	Pg. 8
Getting a blood sample	Pg. 12
Applying blood and reading results	Pg. 24
Reviewing past results and averages	Pg. 32
Control solution testing	Pg. 38
Caring for your system	Pg. 44
Understanding error and other messages	Pg. 50
Detailed information about your system	Pg. 59
Index	Pg. 64

#### Before you begin

Before using this product to test your blood glucose, carefully read this booklet and the inserts that come with the OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips and OneTouch<sup>®</sup> Select<sup>™</sup> and OneTouch<sup>®</sup> Select<sup>™</sup> High Control Solutions. Take note of warnings and cautions throughout this booklet, which are identified with △. Many people find it helpful to practice the test with control solution before testing with blood for the first time. See page 38 in the Control solution testing section.

#### Intended use

The OneTouch<sup>®</sup> Select<sup>™</sup> Blood Glucose Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood. The OneTouch<sup>®</sup> Select<sup>™</sup> System is intended for self-testing outside the body (*in vitro* diagnostic use) by people with diabetes at home and by health care professionals in a clinical setting as an aid to monitor the effectiveness of diabetes control. It should not be used for the diagnosis of diabetes or for testing newborns.

#### **Test principle**

Glucose in the blood sample mixes with special chemicals in the test strip and a small electric current is produced. The strength of this current changes with the amount of glucose in the blood sample. Your meter measures the current, calculates your blood glucose level, displays the result, and stores it in its memory.

#### The OneTouch<sup>®</sup> Select<sup>™</sup> Blood Glucose Monitoring System

Your new OneTouch<sup>®</sup> Select<sup>™</sup> System Kit includes:

- a OneTouch<sup>®</sup> Select<sup>™</sup> Meter (battery included)
- b OneTouch<sup>®</sup> Lancing Device If another type of lancing device was included, see the separate instructions that came with that lancing device.
- c OneTouch<sup>®</sup> AST<sup>™</sup> Clear Cap
- d OneTouch<sup>®</sup> UltraSoft<sup>®</sup> Sterile Lancets
- e Carrying Case

You will also need OneTouch<sup>®</sup> Select<sup>™</sup> Control Solution and OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips.

#### Control solution and test strips are sold separately.

Documents in your kit include this owner's booklet, a starter guide, and a warranty registration card. If any items are missing from your kit, call Customer Service.

▲ WARNING: Keep the meter and testing supplies away from young children. Small items such as the battery door, battery, test strips, lancets, protective disks on the lancets, and control solution vial cap are choking hazards. Do Not ingest or swallow any items.

#### Setting up your system

## Getting to know your OneTouch<sup>®</sup> Select<sup>™</sup> Blood Glucose Meter and test strips

Select or change information

#### Turning your meter on

There are two ways to turn your meter on:

To perform a test, insert a test strip as far as it will go. The meter will briefly perform system checks, then the display will turn on.

or,

With the meter turned off, press and hold off for two seconds to access MAIN MENU.

▲ CAUTION: If you see any light areas within the black start-up screen, there may be a problem with the meter. Call Customer Service.

Check that the screen shows solid black for two seconds. If it does, the display is working properly. If the meter does not power on, try changing the meter battery. See pages 48–49.

#### Turning your meter off

There are several ways to turn your meter off:

- Press and hold 👁 for three seconds.
- Your meter will turn off by itself if left alone for two minutes.
- Go to MAIN MENU and press  $\Delta$  or  $\nabla$  to highlight METER OFF, then press @.
- Before or after completing a test, remove the test strip. If you advance from the test result screen to the MAIN MENU by pressing , removing the test strip will not turn the meter off. Use one of the three methods above.

#### Setting the meter language, date and time

You can change many of the settings that came pre-set with your meter. Before using your meter for the first time or if you change the meter battery, you should check and update these settings. Make sure you complete steps 1 through 6 below to ensure your desired settings are saved.

1. Turn the meter on, see page 2

#### 2. Get to the SETTINGS screen

When using the meter for the first time, or after changing the meter battery, you will automatically start in the SETTINGS screen.

In other cases, from the MAIN MENU, press  $\triangle$  or  $\nabla$  to select SETTINGS. To confirm your selection, press  $\mathfrak{B}$ .

Setting up your system

The pre-set language, date, time and unit of measure appear on the display. If they are correct, press 👁 to return to the MAIN MENU.

If the language, date, or time are not correct, select CHANGE and press P.

The display now shows the LANGUAGE screen.

#### 3. Choose a language

Now press  $\triangle$  or  $\nabla$  to highlight the language of your choice. To confirm your selection, press O.

The display now shows the DATE AND TIME screen.

**NOTE:** If you select the wrong language, press **29** and hold for three seconds to turn the meter off. Then, turn the meter back on and re-start from step 1.

4. Set the date
In the DATE AND TIME screen,
press △ or ▼ to change the year.
To confirm your selection, press ④.

Press  $\triangle$  or  $\nabla$  to change the month. To confirm your selection, press  $\Im$ .

Press  $\triangle$  or  $\nabla$  to change the day. To confirm your selection, press @.

DATE AND TIME
DATE:
TIME: 7 : 56 PM
1 100 1 11
<b>FEB</b> / 27 / 2007
7 : 56 PM
FEB / 27 / 2007
TIME:
7 : 56 PM

5. Set the time of day
Press ▲ or ▼ to set the hour.
To confirm your selection, press ④.

Press  $\triangle$  or  $\nabla$  to set the minutes. To confirm your selection, press @.

Press  $\triangle$  or  $\nabla$  to set AM or PM. To confirm your selection, press O.

The display now shows the SETTINGS screen.

**NOTE:** The date and time formats have been pre-set. You cannot change these settings.

DATE AND TIME		
DATE:		
FFB 1 27	1 2007	
	1 2001	
TIME		
7 · 56 PM		
1 . JU I IVI		

DATEAN	ID TIME
DATE:	
FEB / 27	1 2007
TIME	
7 . 56	DIVI
7 . 30	191

#### 6. Confirm your settings

OK will be highlighted near the bottom of the screen. If your settings are correct, press **1** to confirm and save the settings and return to the MAIN MENU.

▲ WARNING: The unit of measure mg/dL must be displayed here. If your display shows mmol/L rather than mg/dL, contact Customer Service. You cannot change the unit of measure. Use of the wrong unit of measure may cause you to misinterpret your blood glucose level, and may lead to incorrect treatment.

If you want to cancel your settings and start the settings process over again, press  $\triangle$  or  $\nabla$  to highlight CHANGE and press 3. You will be returned to the LANGUAGE screen.

#### Starting the test process

Have these things ready when you test your blood glucose level:

- OneTouch<sup>®</sup> Select<sup>™</sup> Meter
- OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips (sold separately)
- Lancing device
- Sterile lancets with protective disks

#### NOTE:

- Use only OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips with your OneTouch<sup>®</sup> Select<sup>™</sup> Meter.
- Make sure your meter and test strips are about the same temperature before you test.
- Testing must be done within the operating temperature range (50–111°F). For the most accurate results, try to test as close to room temperature (68–77°F) as you can.
- OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips are for single use only. Never re-use a test strip that had either blood or control solution applied to it.

▲ CAUTION: If you cannot test due to a problem with your testing supplies, contact your health care professional or Customer Service. Failure to test could delay treatment decisions and lead to a serious medical condition.

#### **CAUTION:** The

test strip vial contains drying agents that are harmful if inhaled or swallowed and may cause skin or eye irritation.

## **1.** Check the code on the test strip vial before inserting the test strip

Code numbers are used to calibrate your meter with the test strips you are using.

**2. Insert a test strip to turn on the meter** Start with the meter off. If you have turned the meter on to change settings or review past results, turn it off. Remove a test strip from its vial and recap the vial immediately. With clean, dry hands, you may touch the test strip anywhere on its surface. **Do Not** bend, cut or modify the test strips in any way. Use each test strip immediately after removing it from the vial.

Insert the test strip into the test port as shown. Make sure the three contact bars are facing you and the arrow is pointing down. Push the test strip in as far as it will go.

After the black start-up screen appears, the meter will display the code from your last test. If a flashing "——" appears instead of a code number, such as when you are first using the meter, follow the instructions in Step 3 to change to a numerical code.

## **3.** Match the code on the meter with the code on the test strip vial

If the code on the meter does not match the code on the test strip vial, press  $\triangle$  or  $\nabla$  to match the code number on the test strip vial. The new code number will flash on the display for three seconds, then briefly stop flashing, after which the display will advance to the APPLY SAMPLE screen.

If the codes already match, press 👁 to go to the APPLY SAMPLE screen. When you do not make a change after five seconds, the display will advance to the APPLY SAMPLE screen.

The meter is now ready to perform a blood glucose test.

▲ CAUTION: Matching the code on the meter and the code on the test strip vial is essential to obtaining accurate results. Each time you test, check to make sure the code numbers match.

#### NOTE:

- If the APPLY SAMPLE screen appears before you are sure the codes match, remove the test strip, and re-start from step 1, see page 9.
- If you change APPLY SAMPLE to APPLY CONTROL by mistake, press ▼ to change it back to APPLY SAMPLE.

#### Getting a blood sample

#### **Overview of the OneTouch® Lancing Device**

\*The blue cap and OneTouch<sup>®</sup> AST<sup>™</sup> Clear Cap are also used for depth adjustment.

#### **CAUTION:** To reduce the chance of infection:

- Make sure to wash the puncture site with soap and water before sampling.
- Never share a lancet or a lancing device with anyone.
- Always use a new, sterile lancet—lancets are for single use only.
- Keep your meter and lancing device clean. See pages 45–46.

**NOTE:** If you do not have a OneTouch<sup>®</sup> Lancing Device, please refer to the instructions that came with your lancing device.

#### Choosing the right sampling site at the right time

The OneTouch<sup>®</sup> Select<sup>™</sup> Meter allows you to sample blood from your fingertip, forearm or palm. Forearm and palm sampling is also referred to as "alternate site testing" (AST). At times, results obtained at the forearm or palm may be different from a fingertip measurement. Talk to your health care professional before you begin using your forearm or palm for sampling.

If you are testing:	Use blood sample from your:
Routinely before meals	Fingertip, forearm, or palm
<ul> <li>Prior to or more than two hours after:</li> <li>a meal</li> <li>a rapid-acting insulin injection or insulin pump bolus</li> <li>exercise</li> </ul>	Fingertip, forearm, or palm
<ul> <li>When your blood glucose is changing rapidly, such as:</li> <li>within two hours after a meal</li> <li>within two hours after a rapid-acting insulin injection or insulin pump bolus, or</li> <li>during or within two hours after exercise</li> </ul>	Fingertip
When you are concerned about the possibility of hypoglycemia (low blood sugar)	Fingertip

#### **CAUTION: Do Not** test on your forearm or palm when:

- You think your blood glucose is rapidly falling, such as within two hours of exercise or a rapid-acting insulin injection or insulin pump bolus. Testing with a fingertip sample may identify hypoglycemia or an insulin reaction sooner than testing with a forearm or palm sample.
- It has been less than two hours after a meal, a rapid-acting insulin injection or insulin pump bolus, physical exercise, or you think your glucose level is changing rapidly.
- You are concerned about the possibility of hypoglycemia or an insulin reaction, such as when driving a car. This is especially important if you suffer from hypoglycemia unawareness (lack of symptoms to indicate an insulin reaction).

Remember: Consult with your health care professional before using your forearm or palm for testing.

Choose a different puncture site each time you test. Repeated punctures in the same spot may cause soreness and calluses.

If bruising occurs at an alternate site or you have difficulty getting a sample, consider sampling from a fingertip instead. You may want to review the choice of sites with your health care professional.

#### Preparing your sample site

Before you test your blood glucose, wash your hands and forearm (if applicable) thoroughly with warm, soapy water. Rinse and dry.

Lancing and sampling from your fingertip

1. Remove the blue cap by snapping it off

2. Insert a sterile lancet into the OneTouch® Lancing Device Insert the lancet into the holder and push in firmly. Twist the protective disk until it separates from the lancet and save the disk for later use. **Do Not** twist the lancet.

3. Replace the blue cap by snapping it back on

#### 4. Adjust the depth setting

The OneTouch<sup>®</sup> Lancing Device has nine puncture depth settings, numbered 1 through 9. The smaller numbers are for a shallower puncture, and the larger numbers are for a deeper puncture. Shallower punctures work for children and most

adults. Deeper punctures work well for people with thick or callused skin. Twist the blue cap until the correct setting appears.

**NOTE:** A shallower puncture may be less painful. Try a shallower setting first and increase the depth until you find the one deep enough to get a sufficient blood drop for testing.

#### 5. Cock the OneTouch® Lancing Device

Slide the cocking control back until it clicks. If it does not click, it may have been cocked when you inserted the lancet.

#### 6. Puncture your finger

Hold the OneTouch<sup>®</sup> Lancing Device firmly against the side of your finger. Press the release button. Remove the OneTouch<sup>®</sup> Lancing Device from your finger.

#### 7. Get a round drop of blood

Gently squeeze and/or massage your fingertip until a round drop of blood forms on your fingertip.

If the blood smears or runs, **Do Not** use that sample. Wipe the area and gently squeeze another drop of blood or puncture a new site.

#### Lancing and sampling from an alternate site

Sampling from your palm or forearm allows you to use your fingertips less often. You may find that obtaining a blood sample from an alternate site is less painful than using a fingertip. Getting a blood sample from your forearm or palm is different than getting a sample from your fingertips.

#### Forearm sampling

Choose a fleshy area of the forearm away from bone, visible veins and hair. Sometimes there is less blood flow to the forearm than to the fingertips. To help you get a large enough drop of blood, you may gently massage or apply warmth to the site to increase blood flow.

#### Palm sampling

Choose a fleshy area on the palm below your thumb or pinky finger. Select a spot with no visible veins and away from deep lines which may cause your blood sample to smear.

Forearm

The OneTouch<sup>®</sup> AST<sup>™</sup> Clear Cap is used for forearm and palm sampling only. Replace the blue cap with the OneTouch<sup>®</sup> AST<sup>™</sup> Clear Cap.

1. Remove the blue cap by snapping it off

## 2. Insert a sterile lancet into the OneTouch® Lancing Device

Insert the lancet into the holder and push in firmly. Twist the protective disk until it separates from the lancet and save the disk for later use. **Do Not** twist the lancet.

## **3.** Install the OneTouch<sup>®</sup> AST<sup>™</sup> Clear Cap on the OneTouch<sup>®</sup> Lancing Device by snapping it on

#### 4. Adjust the depth setting

You may have to adjust the OneTouch<sup>®</sup> Lancing Device to a deeper setting to get a large enough drop of blood from your forearm or palm. Twist the OneTouch<sup>®</sup> AST<sup>™</sup> Clear Cap toward the larger numbers to increase the depth.

#### 5. Cock the OneTouch<sup>®</sup> Lancing Device

Slide the cocking control back until it clicks. If it does not click, it may have been cocked when you inserted the lancet.

6. Puncture your forearm or palm Firmly press and hold the lancing device against your forearm or palm for a few seconds. Wait until the skin surface under the OneTouch® AST<sup>™</sup> Clear Cap changes color (as blood collects beneath the skin). This tells you there is enough blood flow for a good sample. Then press the release button while continuing to apply pressure. Keep holding the lancing device against your skin until a round drop of blood forms under the cap.

**7.** Remove the OneTouch® Lancing Device Carefully lift the lancing device away from your skin. **Do Not** smear the blood sample.

Forearm

# Getting a blood sample

#### NOTE:

- You may need to wait a little longer to get a large enough drop of blood from the forearm or palm. **Do Not** squeeze the site excessively.
- If the sample drop of blood runs or spreads due to contact with hair or with a line in your palm, **Do Not** use that sample. Try puncturing again in a smoother area.
- Remember: You may have to adjust the lancing device to a deeper setting to get a large enough drop of blood.

#### Applying blood and reading results

Once you have a blood sample and your meter shows the APPLY SAMPLE screen, you are ready to obtain a blood glucose result. If your meter does not show the APPLY SAMPLE screen, remove the unused test strip and re-start the test process. See page 8.

#### 1. Prepare to apply the sample

Keeping your finger extended and steady, move the meter and test strip toward the blood drop.

Fingertip

**Do Not** apply blood on the top of the test strip.

**Do Not** hold the meter and test strip underneath the blood drop. This may cause blood to run into the test port and damage the meter.

When applying a drop of blood from your forearm or palm, keep your palm or forearm steady and bring the top edge of the test strip to the drop of blood with your other hand.

Forearm

#### 2. Apply the sample

Line up the test strip with the blood drop so that the narrow channel on the edge of the test strip is almost touching the edge of the blood drop.

Gently touch the channel to the edge of the blood drop.

Be careful not to push the test strip against your fingertip or the test strip may not fill completely.

#### NOTE:

- **Do Not** smear or scrape the drop of blood with the test strip.
- **Do Not** apply more blood to the test strip after you have moved the drop of blood away.
- Do Not move the test strip in the meter during a test.

**CAUTION:** You may get an ERROR 5 message or an inaccurate result if the blood sample does not fill the confirmation window completely. See page 56. Discard the test strip and re-start the test process.

#### **3.** Wait for the confirmation window to fill completely

The blood drop will be drawn into the narrow channel and the confirmation window should fill completely.

When the confirmation window is full, this means you have applied enough blood. Now you can move the test strip away from the blood drop and wait for the meter to count down from 5 to 1.

Full

#### Not Full

#### 4. Read your result on the meter

Your blood glucose level appears on the display, along with the unit of measure, and the date and time of the test. Blood glucose results are automatically stored in the meter's memory.

#### (Example)

**WARNING:** If mg/dL does not appear with the test result, call Customer Service. Use of the wrong unit of measure may cause you to misinterpret your blood glucose level, and may lead to incorrect treatment.

**CAUTION:** If you test at the low end of the operating range (50°F) and your glucose is high (over 180 mg/dL), the reading on your meter may be lower than your actual glucose. In this situation, repeat the test in a warmer environment with a new test strip as soon as possible.

#### **Error messages**

If you get an ERROR message on your screen rather than a result, see pages 50-58.

#### **Unexpected test results**

Refer to these cautions  $\triangle$  whenever your test results are lower than, higher than, or not what you expect.

#### **CAUTION:** Dehydration and low glucose results

You may get false low glucose results if you are severely dehydrated. If you think you are severely dehydrated, contact your health care professional immediately.

#### ▲ CAUTION: Low glucose results

If your test result is lower than 70 mg/dL or is shown as LOW GLUCOSE, it may mean hypoglycemia (low blood glucose). This may require immediate treatment according to your health care professional's recommendations. Although this result could be due to a test error, it is safer to treat first, then do another test.

#### **CAUTION:** High glucose results

If your test result is higher than 180 mg/dL, it may mean hyperglycemia (high blood glucose). If you are uncertain about this test result, consider re-testing. Your health care professional can work with you to decide what actions, if any, you should take if your results are higher than 180 mg/dL.

If your meter displays HIGH GLUCOSE, you may have a very high blood glucose level (severe hyperglycemia) exceeding 600 mg/dL. Re-check your glucose level. If the result is HIGH GLUCOSE again, this may indicate a severe problem with your blood glucose control and it is important you obtain and follow instructions from your health care professional without delay.

#### **CAUTION:** Repeated unexpected glucose results

- If you continue to get unexpected results, check your system with control solution. See Control solution testing, pages 38–43.
- If you are experiencing symptoms that are not consistent with your blood glucose results and you have followed all instructions in this booklet, call your health care professional. Never ignore symptoms or make significant changes to your diabetes control program without speaking to your health care professional.

#### ▲ CAUTION: Unusual red blood cell count

A hematocrit (percentage of your blood that is red blood cells) that is either very high (above 55%) or very low (below 30%) can cause false results.
# After getting a result

Once you have read your result, you may:

- Review your meter memory by pressing 👁 to go to MAIN MENU, see page 32, or
- Turn the meter off by removing the test strip.

#### Removing the used lancet

Remove the lancing device cap by snapping it off. **Cover the exposed lancet tip before removing the lancet.** Place the lancet protective disk on a hard surface. Push the lancet tip into the disk. Remove the lancet and place it in a container for sharp objects. Replace the cap.

#### Disposing of the used lancet and test strip

It is important to discard the used lancet carefully after each use to avoid unintended lancet stick injuries. Used test strips and lancets may be considered biohazardous waste in your area. Be sure to follow your local regulations for proper disposal.

# **Reviewing past results and averages**

If you have just completed a test, press I to get to the MAIN MENU screen. If your meter is off, press and hold I to turn it on. From the MAIN MENU screen you can choose:

- LAST RESULT to view your most recent result,
- ALL RESULTS to review up to 350 of your most recent results four at a time, or

• AVERAGES to view the average of your results.

Press  $\Delta$  or  $\nabla$  to highlight LAST RESULT, ALL RESULTS, or AVERAGES.

To confirm your selection, press 👁.

# Last result

The meter will display your most recent result. This result will include the unit of measure and the date and time of the test. CONTROL TEST will also appear if this is a control solution test. Press 👁 to return to MAIN MENU.

# All results

The meter will display four results at a time, in the order the tests were done, starting with the most recent. For each test result, the meter will display the date and time of the test. Results may also contain the following symbols:

- ${\rm HI}\,$  if the result was above 600 mg/dL
- ${\rm LO}$  if the result was below 20 mg/dL
- C if the result is from a control solution test

The meter stores a maximum of 350 blood glucose or control solution test results. When the memory is full, the oldest result is dropped as the newest is added.

To view details of an individual result, press  $\triangle$  or  $\nabla$  to highlight the result you want, then press 0. To return to the list of all results from an individual result, press 0.

To view more recent results, continue to press  $\Delta$  after the top result on the display is highlighted.

To view older results, press  $\nabla$  after MENU is highlighted. Pressing and holding  $\Delta$  or  $\nabla$  allows you to move more quickly through the results. To return to the main menu, highlight MENU, then press  $\mathfrak{B}$ .

#### Averages

The meter will display your averages for each of the 7-, 14-, and 30-day periods leading up to the current date.

In result averages, a HIGH GLUCOSE result is counted as 600 mg/dL, and a LOW GLUCOSE result as 20 mg/dL. Control solution results are not part of your averages.

To return to the main menu, press 👁.

## NOTE:

- Result averages provide information from past results. **Do Not** use result averages to make immediate treatment decisions.
- If you change your date setting, your averages may change too. The meter calculates averages based on the 7-, 14-, and 30-day periods ending on the current date setting.

To see averages for different parts of the day or over a different number of days, you can use OneTouch<sup>®</sup> Diabetes Management Software and your home computer. See the next section, on page 36.

#### Downloading results to a computer

You can use your meter with OneTouch® Diabetes Management Software (DMS) for storing your records and to help you spot patterns for planning meals, exercise, and medication. OneTouch® DMS puts information downloaded from the meter into charts and graphs.

#### 1. Obtain the required software and cable

For order information and to learn more about OneTouch<sup>®</sup> Diabetes Management Software, visit www.OneTouchDiabetesSoftware.com.

#### 2. Install the software on a computer

Follow the installation instructions provided with  $\mathsf{OneTouch}^{\circledast}\,\mathsf{DMS}$  and the  $\mathsf{OneTouch}^{\circledast}$  Interface Cable.

**WARNING:** To avoid a possible shock, Do Not insert a test strip when the meter is connected to a computer with the OneTouch<sup>®</sup> Interface Cable.

## 3. Get ready to transfer readings

Make sure that the meter is turned off. Connect the OneTouch<sup>®</sup> Interface Cable to the COM or USB port on your computer. Then connect the other end of the cable to the meter data port.

#### 4. Transfer data

Follow the instructions provided with  $\mathsf{OneTouch}^{\circledast}\mathsf{DMS}$  to download results from the meter.

Once the command to start the download is sent from the computer to the meter, the meter display will show "PC" indicating that the meter is in communication mode. You will not be able to perform a test when the meter is in communication mode.

# **Control solution testing**

You may check your meter and test strips using either the OneTouch<sup>®</sup> Select<sup>™</sup> or OneTouch<sup>®</sup> Select<sup>™</sup> High Control Solutions.

OneTouch<sup>®</sup> Select<sup>™</sup> and OneTouch<sup>®</sup> Select<sup>™</sup> High Control Solutions contain a known amount of glucose and are used to check that the meter and the test strips are working properly. Remember which control solution you used.

Do a control solution test:

- to practice the test process instead of using blood,
- once a week,
- whenever you open a new vial of test strips,
- if you suspect the meter or test strips are not working properly,
- if you have had repeated unexpected blood glucose results (as described on page 28), or
- if you drop or damage the meter.

**CAUTION: Do Not** swallow control solution; it is not for human consumption.

**Do Not** apply control solution to the skin or eyes as it may cause irritation.

# NOTE:

- Use only OneTouch<sup>®</sup> Select<sup>™</sup> or OneTouch<sup>®</sup> Select<sup>™</sup> High Control Solution with your OneTouch<sup>®</sup> Select<sup>™</sup> Meter.
- Control solution tests must be done at room temperature (68–77°F). Make sure your meter, test strips, and control solution are at room temperature before testing.

## Performing a control solution test

Start with the meter off. If you have turned the meter on to change settings or review past results, turn it off.

**1.** Check the code on the test strip vial before inserting the test strip

2. Insert a test strip to turn on the meter Make sure the three contact bars are facing you and the arrow is pointing down. Push the test strip in as far as it will go. **Do Not** bend the test strip.

# **3.** Match the code on the meter with the code on the test strip vial

If the code on the meter does not match the code on the test strip vial, press △ or ♥ to match the code number on the test strip vial. The new code number will flash on the display for three seconds, then briefly stop flashing, after which the display will advance to the APPLY SAMPLE screen.

If the codes already match, press to go to the APPLY SAMPLE screen. When you do not make a change after five seconds, the display will advance to the APPLY SAMPLE screen.

4. Mark the test as a control solution test Press ▲ to change APPLY SAMPLE to APPLY CONTROL. You must mark the test before you apply control solution. Once you have completed the test, you cannot change the marking.

The meter is now ready to perform a control solution test.

## 5. Prepare and apply control solution

Shake the control solution vial before each test. Remove the cap and squeeze the vial to discard the first drop. Then wipe the tip with a clean tissue or cloth. Hold the vial upside down and gently squeeze out a hanging drop.

Touch and hold the hanging drop of control solution to the narrow channel in the **top edge** of the test strip. Make sure the confirmation window fills completely. Control solution should not be applied to the flat face of the test strip.

## 6. Read your result

When the confirmation window is full, the meter will count down from 5 to 1.

Your control solution result will then appear on the display, along with the date, time, unit of measure, and the words CONTROL TEST. The control solution results can be viewed in the list of past results, but are not counted in your result averages.

# **CONFIRMATION WINDOW**

## 7. Check if the result is in range

Each vial of test strips has both OneTouch<sup>®</sup> Select<sup>™</sup> and OneTouch<sup>®</sup> Select<sup>™</sup> High Control Solution ranges printed on its label. Compare the result displayed on the meter to either the OneTouch<sup>®</sup> Select<sup>™</sup> or OneTouch<sup>®</sup> Select<sup>™</sup> High Control Solution range printed **on the test strip vial**, depending on the type of control solution you used. If the results you get are not within this range, the meter and test strips may not be working properly. Repeat the control solution test.

Out-of-range results may be due to:

- not following the instructions detailed on pages 38-42,
- expired or contaminated control solution,
- expired or damaged test strip,
- use of a test strip or control solution past its discard date, or
- a problem with the meter.

(Example ranges)

**NOTE:** It is important to compare your control solution result with the correct control range printed on the test strip vial or your result may appear to be out of range.

▲ CAUTION: The control solution ranges printed on the test strip vial are for OneTouch<sup>®</sup> Select<sup>™</sup> and OneTouch<sup>®</sup> Select<sup>™</sup> High Control Solutions only. They are not recommended ranges for your blood glucose level.

**CAUTION:** If you continue to get control solution test results that fall outside the ranges printed on the test strip vial, **Do Not** use the meter, the test strips, or the control solution. Call Customer Service.

# **Caring for your system**

Your OneTouch<sup>®</sup> Select<sup>™</sup> Blood Glucose Monitoring System does not need any special maintenance.

#### Storing your system

Store your meter, test strips, control solution and other items in your carrying case after each use. Store each item in a cool, dry place below 86°F, but **Do Not** refrigerate. Keep all items away from direct sunlight and heat.

Tightly close the cap on the test strip vial and/or control solution vial immediately after use to avoid contamination or damage. Store test strips only in their original vial. **Do Not** open the test strip vial until you are ready to test.

# Checking for expiration or damage to test strips and control solution

Test strips and control solution have expiration dates printed on their vials. When you first open a test strip or control solution vial, you must record the discard date (date opened plus three months) in the space provided on the label.

**CAUTION: Do Not** use the test strips or control solution after the expiration date printed on the vial or the discard date, whichever comes first, or your results may be inaccurate.

▲ CAUTION: Do Not use your test strips if your vial is damaged or left open to air. This could lead to error messages or tests that read higher than the actual value. Call Customer Service immediately if the test strip vial you just purchased is damaged.

#### **Cleaning your meter**

To clean your meter, wipe the outside with a soft cloth dampened with water and mild detergent. **Do Not** use alcohol or another solvent to clean your meter.

**Do Not** get any liquids, dirt, dust, blood, or control solution inside the meter through the test port or the data port. Never spray cleaning solution on the meter or immerse it in any liquid.

# Cleaning your OneTouch<sup>®</sup> Lancing Device and OneTouch<sup>®</sup> AST<sup>™</sup> Clear Cap

To clean these items, wipe them with a soft cloth dampened with water and mild detergent. **Do Not** immerse the OneTouch<sup>®</sup> Lancing Device in any liquid.

To disinfect these items, prepare a solution of one part household bleach to ten parts water. Wipe the OneTouch<sup>®</sup> Lancing Device with a soft cloth dampened with this solution. Immerse the **caps only** in this solution for 30 minutes. After disinfecting, rinse briefly with water and allow both to air dry.

#### Battery

Your OneTouch<sup>®</sup> Select<sup>™</sup> Meter uses one 3.0 Volt CR 2032 lithium battery (or equivalent). Replacement batteries can be found in most stores where batteries are sold. Your meter comes with the battery already installed.

#### Low meter battery

The meter shows a battery icon ( ) in the lower left corner of the display, or a low battery message to indicate the condition of the meter battery. When the battery icon or the LOW BATTERY/REPLACE BATTERY SOON screen message first appears, there is enough power for a minimum of 100 more tests. You should replace the meter battery as soon as possible.

When your meter displays the LOW BATTERY/ REPLACE BATTERY NOW! screen message, there is not enough battery power remaining to perform a test. You must install a new battery before using your meter.

## **Replacing the battery**

## **1. Remove the old battery** Start with the meter off. Open the battery door and pull up on the battery ribbon.

**WARNING:** To avoid a possible shock, Do Not change the battery while the meter is connected to a computer with the OneTouch<sup>®</sup> Interface Cable.

#### 2. Insert the new battery

With the "+" side facing up toward you, place the battery in the compartment within the fold of the ribbon. Push the battery until it snaps into the battery clasp. Insert the two battery door tabs into the matching holes on the meter, and push down until you hear the door click into place.

If the meter does not power on after you have replaced the meter battery, check that the battery is correctly installed with the "+" side up. If the meter still does not power on, call Customer Service.

#### 3. Check your meter settings

Removing the meter battery will not affect your stored results. However, you may need to re-set your meter settings. See pages 3–7.

# **4.** Dispose of the battery according to your local environmental regulations

Lithium ion batteries contain perchlorate material. Special handling may apply, see California Code of Regulations Title 22, Section 67384.4.

# Understanding error and other messages

The OneTouch<sup>®</sup> Select<sup>™</sup> Meter displays messages when there are problems with the test strip, with the meter, or when your blood glucose levels are higher than 600 mg/dL or lower than 20 mg/dL. Messages do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message.

Message	What it means	What to do	
WARNING	You may have a very	This may require immediate	
LOW GLUCOSE	low blood glucose level (severe hypoglycemia),	health care professional's	
BELOW 20 mg/dL	lower than 20 mg/dL.	recommendations.	
		be due to a test error, it is	
		safer to treat first and then do another test.	

Message	What it means	What to do	
WARNING HIGH GLUCOSE ABOVE 600 mg/dL	You may have a very high blood glucose level (severe hyperglycemia), over 600 mg/dL.	Re-check your glucose level. If the result is HIGH GLUCOSE again, obtain and follow instructions from your health care professional without delay.	

WADNING	Your control solution	Test again with control solution
- WARINING	test result is below	and a new test strip. If the result
LOW GLUCOSE	20 mg/dL.	is LOW GLUCOSE again, <b>Do Not</b>
BELOW/20 mg/dL		use the meter. Contact Customer
		Service.
CONTROL TEST		nd

Message	What it means	What to do
WARNING HIGH GLUCOSE ABOVE 600 mg/dL	Your control solution test result is above 600 mg/dL.	Test again with control solution and a new test strip. If the result is HIGH GLUCOSE again, <b>Do Not</b> use the meter. Contact Customer
CONTROL TEST		Service.

Understanding error and other messages

Message	What it means	What to do	
WARNING TEMPERATURE ERROR TEMPERATURE TOO LOW SEE OWNER'S BOOKLET	Meter is too cold (below 50°F) to work correctly.	Wait a few minutes and insert a new test strip. If you do not get another TEMPERATURE ERROR message, the meter is now within the operating range.	
WARNING ERROR 1 METER PROBLEM CALL CUSTOMER SERVICE	There is a problem with the meter.	<b>Do Not</b> use the meter. Contact Customer Service.	and
WARNING ERROR 2 METER OR STRIP PROBLEM RETEST WITH A NEW STRIP	Error message could be caused either by a used test strip or a problem with the meter.	Repeat the test with a new test strip; see pages 24–27. If this message continues to appear, contact Customer Service.	other messages

Message	What it means	What to do
WARNING	The sample was	Repeat the test with a new test
ERROR 3	applied before the	strip. Apply a blood or control
	meter was ready.	Solution sample only after
RETEST WITH A		appears on the display. If this
NEW STRIP		message continues to appear,
		contact Customer Service.

	WARNING ERROR 4	One of the following may apply:	
or es	STRIP PROBLEM	You may have high glucose and have	If you tested in a cool environment, repeat the test in a
Understanding erro and other message	SEE OWNER'S BOOKLET	tested in an environment near the low end of the system's operating temperature range (50–111°F). <i>or</i> ,	warmer environment with a new test strip; see pages 24–27. If the error message appears again, contact Customer Service.

Message	What it means	What to do
WARNING ERROR 4 STRIP PROBLEM SEE OWNER'S BOOKLET	There may be a problem with the test strip. For example, it may have been damaged or moved during testing. <i>or</i> ,	If you tested in a normal or warm environment, repeat the test with a new test strip; see pages 24–27. If the error message appears again, contact Customer Service.
	The sample was improperly applied. <i>or,</i>	If you applied the sample incorrectly, review pages on blood application (see pages 24–27) or control solution testing (see pages 38–43) and repeat the test with a new test strip. If the error message appears again, contact Customer Service.
	There may be a problem with the meter.	If the error message appears again, contact Customer Service.

	Message	What it means	What to do
	WARNING ERROR 5 STRIP PROBLEM OR SAMPLE TOO SMALL RETEST WITH A NEW STRIP	The meter has detected a problem with the test strip. Possible causes are test strip damage or an incompletely filled confirmation window.	Repeat the test with a new test strip. Refer to information on blood application (see pages 24–27) or control solution testing (see pages 38–43).
Understanding error and other messages	WARNING D LOW BATTERY REPLACE BATTERY SOON SEE OWNER'S BOOKLET	Meter battery is low but still has enough power to perform a test.	This message will appear when you first turn the meter on and there is enough power for a minimum of 100 more tests. Press 👁 to continue, but replace the battery as soon as possible.

Message	What it means	What to do	
MAY/16/07 2:19 PM 123 mg/dL	Meter battery is low but still has enough power to perform a test.	Once the LOW BATTERY/REPLACE BATTERY SOON message is displayed, this icon will continue to appear on screens until you replace the battery. When the battery icon first appears, there is enough power for a minimum of 100 more tests. Test results will still be accurate, but replace the battery as soon as possible.	
WARNING D LOW BATTERY REPLACE BATTERY NOW! SEE OWNER'S BOOKLET	Meter battery does not have enough power to perform a test.	Replace meter battery.	מוות טרוופו ו

Message	What it means	What to do
	No result in memory, such as the first time use of the meter or after a download of all data to a computer. <i>or</i> , Your meter was unable to recall this result. This result will not be included in result averages.	<ul> <li>You can still perform a blood glucose test and get an accurate result. Contact Customer Service to report this occurrence if this is <b>not</b>:</li> <li>1. Your first time use of the meter, or</li> <li>2. After you just downloaded results from your meter.</li> </ul>
ALL RESULTS	Your meter was unable	You can still perform a blood

Understanding error and other messages

Your meter was unabl to recall this result. This result will not be included in result averages. You can still perform a blood glucose test and get an accurate result, but contact Customer Service to report this occurrence.

# Detailed information about your system

#### Comparing meter and laboratory results

Test results with the OneTouch<sup>®</sup> Select<sup>™</sup> Meter are plasma-calibrated. This helps you and your health care professional to compare your meter results with laboratory tests. If you have been using another type of meter—one that provides whole blood-calibrated results—you may notice that your test results with the OneTouch<sup>®</sup> Select<sup>™</sup> Meter are approximately 12% higher.

OneTouch<sup>®</sup> Select<sup>™</sup> Meter test results and laboratory test results both are expressed in plasma-equivalent units. However, your meter result may differ from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way.

Your OneTouch<sup>®</sup> Select<sup>™</sup> Meter glucose value is considered accurate when it is within ±20% of the laboratory measurement. There are some specific situations that could cause a difference of more than ±20%:

• You have eaten recently. The blood glucose level from blood obtained from a fingertip can be up to 70 mg/dL higher than blood drawn from a vein (venous sample) used for a lab test.<sup>1</sup>

<sup>1</sup> Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), *Tietz Textbook of Clinical Chemistry*. Philadelphia: W.B. Saunders Company (1994), 959.

- Your hematocrit (percentage of your blood that is red blood cells) is high (above 55%) or low (below 30%).
- You are severely dehydrated.
- You tested at a temperature near the low end of the operating range (50°F) and you get a high glucose result (i.e., greater than 180 mg/dL). In this situation, repeat the test in a warmer environment with a new test strip as soon as possible.

For accuracy and precision data and for important information on limitations, see the insert that comes with your test strips.

To maximize your chances of an accurate comparison between meter and laboratory results, follow a few basic guidelines:

# Before going to the lab

- Perform a control solution test to make sure the meter is working properly.
- **Do Not** eat for at least eight hours before you test your blood.
- Take your meter with you to the lab.

# While at the lab

- Conduct your meter test within 15 minutes of the lab test.
- Use only fresh, capillary blood obtained from the fingertip.
- Follow all instructions in this owner's booklet for performing a blood glucose test with your meter.

Detailed information about your system

# **Technical specifications**

Reported result range	20-600 mg/dL
Calibration	Plasma-equivalent
Sample	Fresh capillary whole blood
Test time	5 seconds
Assay method	Glucose oxidase biosensor
Meter power source	One replaceable 3.0 Volt CR 2032 lithium battery (or equivalent)
Unit of measure	mg/dL
Memory	350 blood glucose or control solution test results
Automatic shutoff	2 minutes after last action
Size	3.54 x 2.19 x 0.85 inches
Weight	Approximately 1.8 ounces, with battery

Detailed information about your system

Operating ranges	Temperature: 50–111°F Relative humidity: 10–90% Altitude: up to 10,000 feet Hematocrit: 30–55%
Battery ratings	$1 \times 3.0 \text{ V}$ d.c., 60 mA (1 x CR 2032 battery) direct current

#### Symbols

Cautions and Warnings: Refer to safety-related notes in the owner's booklet and inserts that came with your meter and testing supplies.

Low	battery

Direct current

# Guarantee

LifeScan guarantees that the OneTouch<sup>®</sup> Select<sup>™</sup> Meter will be free of defects in material and workmanship for three years, valid from the date of purchase. The guarantee extends only to the original purchaser and is not transferable.

## **Electrical and safety standards**

This meter complies with CISPR 11: 2003, Class B (Radiated Only). Emissions of the energy used are low and not likely to cause interference in nearby electronic equipment.

The meter has been tested for immunity to Level 3 electrostatic discharge as specified in IEC 61000-4-2.

This meter has been tested for immunity to radio frequency interference over the frequency range 80MHz to 2.5GHz at 3V/m as specified in IEC 61000-4-3.

Degree of protection rating: IP31

CAN/CSA C22.2 61010-1:04, UL 61010-1:04, IEC 61010-1 and LEC 61010-2-101

# Index

ALL RESULTS AVG
Alternate site testing procedure19
Averages
Battery 1, 47
Battery, icon on display 47,56, 57
Black start-up screen2, 9
Buttons on meter 1
Cable, OneTouch <sup>®</sup> Interface
Cleaning your meter and lancing device 45, 46
Code9, 10, 39
Confirmation window 1, 27, 41
Connecting to a computer
Control solutions, result ranges
Control solutions, storing44
Control solutions, testing
Control solutions, types
Data port 1, 37
Date setting 5
Dehydration
Discard date, control solutions45
Discard date, test strips45
Display check 2, 9
Downloading results to a computer 36, 37
Error messages

Expiration date, control solutions
Expiration date, test strips
Fingertip testing procedure
Forearm
Forearm or palm sampling.
important information 13 14 23
Guarantee 62
HIGH GUICOSE Control Test message 52
Hunorghucomia
Hypoglycemia
Infection, reduce the chance12
Intended use
Kit componentsiv
Lab results, comparing with meter results59
Lancet, disposal
Lancet, inserting in lancing device 16, 20
Lancet, protective disk 12, 16, 20, 31
Lancing 15-23, 46
Lancing device and capiv, 12, 16, 46
Lancing device,
OneTouch <sup>®</sup> AST <sup>™</sup> Clear Cap iv, 12, 21, 46
Lancing device, depth setting 12, 16, 21
Language, choosing in set up4

LOW BATTERY message
LOW GLUCOSE Control Test message51
LOW GLUCOSE message 29, 34, 50
MAIN MENU screen
mg/dL7, 28, 61
mmol/L
Palm 13-14, 19-22
PC communication mode
Plasma calibration
Repeated unexpected glucose results 30
Result averages, 7-, 14-, and 30-day34
Result screen
Result, last
Results, all
Results, reviewing past
Settings
Software, OneTouch® Diabetes Management . 36
Start-up screen
Storing your control solutions
Storing your meter
Storing your test strips 44
Symbols 1, 33, 62
Technical specifications
Temperature
TEMPERATURE ERROR messages 52, 53
Testing procedure8

Test strip, applying drop of blood 24, 25
Test strip, confirmation window 1, 27, 41
Test strip, contact bars 1, 9, 39
Test strip, disposal
Time setting6
Turning meter off2
Turning meter on for reviewing results 2, 32
Turning meter on for testing
Unexpected test results
Unit of measure

Your target ranges
from your health care professional

### **Patent information**

Density dens

The system described herein is covered by one or more of the following U.S. patents: 5,708,247, 5,951,836, 6,241,862, 6,284,125, and 7,112,265. Use of the monitoring device included herein is protected under one or more of the following U.S. patents: 6,413,410, 6,733,655. Purchase of this device does not act to grant a use license under these patents. Such a license is granted only when the device is used with OneTouch<sup>®</sup> Select<sup>™</sup> Test Strips. No test strip supplier other than LifeScan is authorized to grant such a license. The accuracy of results generated with LifeScan meters using test strips manufactured by anyone other than LifeScan has not been evaluated by LifeScan.
#### ¡Bienvenido a la familia OneTouch®!

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Manufactured for / Fabricado para: LifeScan, Inc. Milpitas, CA 95035 USA © 2007 LifeScan, Inc. Rev. date / Fecha de rev.: 8/2007