

miniland
BABY

little big ideas

thermotalk plus

The thermometer that gives out a voice message and with pulse function

REF: 89068



- ⊙ Ultrafast thermometer with temperature reading in 1 second
- ⊙ It takes the temperature on the baby's forehead and in the baby's ear
- ⊙ Fever indicator
- ⊙ With pulse reading in 15 seconds, to control the parents' heartbeat
- ⊙ It memorizes the last 15 readings with additional information about the date and type of reading (forehead or ear)
- ⊙ Room temperature reading
- ⊙ It gives the temperature with a voice message, in 6 languages, Spanish, English, German, French, Italian and Russian
- ⊙ Silent mode
- ⊙ With clock and calendar

thermotalk plus

The thermometer that gives out a voice message and with pulse function

miniland
BABY

little big ideas



Technical Features

- Technology: Infra red sensor
- Reading time: 1 second
- Batteries: 2 x AAA. Power: DC 3 V
- Measuring range: 32,0 °C ~ 42,9 °C (89,6 °F ~ 109,3 °F)
- Environmental temperature: 5,0 °C ~ 59,9 °C (41,0 °F ~ 139,8 °F)
- Measuring precision:
 - +/- 0,2 °C (or +/- 0,4 °F)
from 35,5 °C ~ 42,0 °C (or 95,9 °F ~ 107,6 °F)
 - +/- 0,3 °C (+/- 0,5 °F)
from 35,0 °C ~ 35,4 °C (or 89,6 °F ~ 95,8 °F)
from 41,2 °C ~ 42,9 °C (or 107,7 °F ~ 109,3 °F)
- Environmental temperature: +/- 2 °C (4 °F)
- Working environment:
 - 10 °C ~ 40 °C (50,0 °F ~ 104,0 °F) for ear
 - 15 °C ~ 40 °C (59,0 °F ~ 104,0 °F) for forehead
 - Relative humidity: 20% RH ~ 85% RH
- Storage environment:
 - 10 °C to 55 °C (14 °F a 131 °F) Relative humidity 20% RH ~ 85% RH
- Precision / Repeated: +/- 0,3 °C (+/- 0,5 °F)
- External dimensions: 127 (length) x 39 (width) x 63 (depth) mm
- Weight: 67 g (without batteries)

Functions

- Use: people/environment
- Voice communication: report with reading results, room temperature and time
- Option 6 different languages: Spanish, English, German, French, Italian and Russian
- Readings: ear/forehead
- N° of memorised readings: 15
- High temperature alert
- Low battery alert
- Pulse reading
- Back-lit LCD screen